

DRYLAB KNOWLEDGE MANAGEMENT DOCUMENT

UPLC Method for Amlodipine: Aquity BEH C18

by R. Kormány, I. Molnár, H.J. Rieger

Summary

Project Description


Exchange an old HPLC method against a better and faster method.

Analytical Target Profile (ATP)


Separate all the drug substance from all the impurities with a minimum resolution of 1.5. The limit of the quantification must be less than 0.05% of the drug substance (reporting limit according ICH Q3A).

Confirmation Run vs Predicted Run

Predicted Run



Confirmation Run



Conclusions

In reference to the analytical target profile (ATP) the selected working point shows a separation off all impurities from each other and from the drug substance with more than resolution of 1.5 (2.6). The separation is good enough to guarantee the quantification limit 0.05%. The new proposed method has a run time minor than 6 min and it is compared to the old method (60) 10 times faster.

Signature:

Name:

Table of Contents

DryLab Knowledge Management Document	2
Project Description	2
Analytical Target Profile (ATP)	2
Risk Management	2
Documentation of the initial experiments	3
3D Experimental Space	11
Selected Working Point (Set Point)	12
Method Robustness Calculation	13
Confirmation Run vs Predicted Run	14
Conclusions	15
References and Comments	15
Annex	16

DRYLAB KNOWLEDGE MANAGEMENT DOCUMENT

UPLC Method for Amlodipine: Acquity BEH C18

by R. Kormány, I. Molnár, H.J. Rieger

Project Description

Exchange an old HPLC method against a better and faster method.

Analytical Target Profile (ATP)

Separate all the drug substance from all the impurities with a minimum resolution of 1.5. The limit of the quantification must be less than 0.05% of the drug substance (reporting limit according ICH Q3A).

Risk Management

The Critical Quality Attribute (CQA) is: sufficient separation of all the impurities according to the ATP.

Based on prior knowledge and experience with similar projects, the parameters: gradient time tG (3-9 [min]), pH (2 - 2.5 - 3), temperature T (15-45 [°C]) and ternary composition of the eluent B are the most probable potentially critical separation parameters.

Initial Conditions

The three parameters gradient time tG (3-9 [min]), temperature T(15-45 [°C]) and pH (2-2.5-3) were optimized. The stationary phase Acquity HSS C18 was selected based on a previous test plan [1].

min	°C	
tG1: 3.00	T1: 15.00	pH1: 2.00
tG2: 9.00	T2: 45.00	pH2: 2.50
		pH3: 3.00


Column Data	
Name: Acquity BEH C18	
Diameter [mm]:	2.10
Length [mm]:	50.00
Particle size [µm]:	1.70
to [min]:	0.21
FlowRate [mL/min]:	0.50
Injection Volumen [µL]:	1.00

Instrument Data	
Name: N/A	
Dwell volume [mL]:	0.125
Extracol volume [mL]:	0.002
Time constant [s]:	0.10
Wavelength [nm]:	230

Eluent Data	
Eluent A :	
1: 100 % Water (H2O)	
Eluent B :	
1: 100 % Acetonitrile (ACN)	

Documentation of the initial experiments


Reference runs: pH (2.00, 2.50, 3.00)




Run	t _R [min]	T[°C]	pH
2	9.00	15.00	2.00
6	9.00	15.00	2.50
10	9.00	15.00	3.00

Chromatograms of the Reference Runs


Run 10



Run 6




Run 2



Peak Tracking of the References Runs

#	Name	SDev	Run 2		Run 6		Run 10	
			tR [min]	Area	tR [min]	Area	tR [min]	Area
1	ImpD	1.36	0.85	261	0.86	258	0.92	252
2	ImpF	1.73	0.96	354	0.98	353	1.04	340
3	Amlodipine	2.00	1.27	366	1.29	366	1.37	351
4	ImpE	1.83	1.62	324	1.65	322	1.74	311
5	ImpB	0.94	2.58	512	2.59	512	2.65	502
6	ImpG	1.16	2.71	570	2.71	570	2.77	556
7	ImpH	1.55	2.89	484	2.89	482	2.92	467
8	ImpA	0.99	5.13	673	5.13	671	5.16	658
	Σ Areas	1.36	18	3543	18	3534	19	3437


tG-T-Plane 1 (pH 2.00)




Run	tG[min]	T[°C]	pH
1	3.00	15.00	2.00
2	9.00	15.00	2.00
3	3.00	45.00	2.00
4	9.00	45.00	2.00

Chromatograms of tG-T-Plane 1


Run 3




Run 4



Run 1




Run 2




Peak Tracking of the tG-T-Plane 1

#	Name	SDev	Run 1		Run 2		Run 3		Run 4	
			tR [min]	Area						
1	ImpD	0.61	0.73	259	0.85	261	0.72	257	0.85	262
2	ImpF	0.57	0.78	350	0.96	354	0.77	349	0.93	353
3	Amlodipine	0.54	0.93	362	1.27	366	0.92	363	1.25	366
4	ImpE	0.46	1.08	321	1.62	324	1.09	322	1.62	325
5	ImpB	1.29	1.51	502	2.58	512	1.51	517	2.57	519
6	ImpG	0.75	1.62	568	2.71	570	1.51	575	2.44	579
7	ImpH	0.19	1.62	483	2.89	484	1.56	484	2.74	485
8	ImpA	1.08	2.51	666	5.13	673	2.29	675	4.61	686
	Σ Areas	0.64	11	3511	18	3543	10	3542	17	3575

Original Resolution Map Plane 1




tG-T-Plane 2 (pH2 2.50)




Run	tG[min]	T[°C]	pH
5	3.00	15.00	2.50
6	9.00	15.00	2.50
7	3.00	45.00	2.50
8	9.00	45.00	2.50

Chromatograms of tG-T-Plane 2


Run 7




Run 8



Run 5




Run 6




Peak Tracking of the tG-T-Plane 2

#	Name	SDev	Run 5		Run 6		Run 7		Run 8	
			tR [min]	Area						
1	ImpD	0.28	0.73	258	0.86	258	0.73	257	0.85	258
2	ImpF	0.35	0.79	350	0.98	353	0.77	349	0.93	350
3	Amlodipine	0.16	0.93	364	1.29	366	0.92	365	1.25	364
4	ImpE	0.29	1.08	322	1.65	322	1.08	322	1.62	320
5	ImpB	0.87	1.50	508	2.59	512	1.50	520	2.57	515
6	ImpG	0.59	1.62	575	2.71	570	1.50	579	2.43	575
7	ImpH	0.60	1.62	487	2.89	482	1.55	488	2.72	481
8	ImpA	0.72	2.50	673	5.13	671	2.28	681	4.59	683
Σ Areas		0.29	11	3539	18	3534	10	3561	17	3547

Original Resolution Map Plane 2




tG-T-Plane 3 (pH3 3.00)




Run	tG[min]	T[°C]	pH
9	3.00	15.00	3.00
10	9.00	15.00	3.00
11	3.00	45.00	3.00
12	9.00	45.00	3.00

Chromatograms of tG-T-Plane 3


Run 11




Run 12



Run 9




Run 10



Peak Tracking of the tG-T-Plane 3


#	Name	SDev	Run 9		Run 10		Run 11		Run 12		
			tR [min]	Area							
1	ImpD	0.59	0.77	254	0.92	252	0.73	250	0.86	253	
2	ImpF	0.41	0.83	343	1.04	340	0.78	340	0.94	340	
3	Amlodipine	0.40	0.98	355	1.37	351	0.93	352	1.26	352	
4	ImpE	0.39	1.13	314	1.74	311	1.10	312	1.64	312	
5	ImpB	0.67	1.53	504	2.65	502	1.51	495	2.57	500	
6	ImpG	7.47	1.65	653	2.77	556	1.51	548	2.43	557	
7	ImpH	7.51	1.63	392	2.92	467	1.54	477	2.69	461	
8	ImpA	0.58	2.52	660	5.16	658	2.29	650	4.60	658	
Σ Areas			0.56	11	3474	19	3437	10	3423	17	3432

Original Resolution Map Plane 3




3D Experimental Space

Max Crit. Resolution: 3.35



Method Operable Design Region (MODR)

Critical Resolution Level: 1.5



Selected Working Point (Set Point)

The working point from the previously constructed Design Space was selected because it had the highest critical resolution (R_s , crit).


Final Conditions

Column Data		Instrument Data		Additional Column Data	
Name:	Acquity BEH C18	Name:	N/A	Observed t_0 [min]:	213.000
Diameter [mm]:	2.10	Dwell volume [mL]:	0.125	at a Flow Rate of [mL/min]:	0.50
Length [mm]:	50.00	Extracol volume [mL]:	0.002	Pore Diameter [nm]:	10.00
Particle size [μm]:	1.70	Time constant [s]:	0.10	A-value:	0.80
t_0 [min]:	213.00	Wavelength [nm]:	230	Plate Number:	0.00
FlowRate [mL/min]:	0.50				
Injection Volumen [μL]:	1.00				


Status		Gradient Table		
8.10	tG [min]	Time [min]	[%B]	Rate [%B/min]
13.50	T [$^{\circ}\text{C}$]	0.00	30.00	
2.10	pH	8.10	90.00	7.41

Pressure [psi]:	6207
Plate Number:	10606 (calculated)
R_s ,crit:	2.96
Crit. Peak Pair:	6, 7
Run Time [min]:	6.00
Eluent Used [mL]:	3.00

Resolution Map for the Selected Working point



Chromatogram Selected Working point



Results Table

#	Name	tR [min]	Area	Avg. k	Width	Resolution
1	ImpD	1.10	259.96	3.22	0.04	3.08
2	ImpF	1.22	352.93	3.31	0.04	7.65
3	Amlodipine	1.51	365.74	3.64	0.04	8.12
4	ImpE	1.82	323.13	3.95	0.04	20.33
5	ImpB	2.69	509.08	5.63	0.05	3.76
6	ImpG	2.87	568.10	6.91	0.05	1.74
7	ImpH	2.96	485.59	5.83	0.05	38.00
8	ImpA	4.90	670.50	8.04	0.06	0.00

Method Robustness Calculation

1.-Frequency distribution:

Out of the 729 model simulation no result was less than 1.52. This is good to guarantee an adequate robustness.

2.-Regression Coefficients:


From the selected set of potentially critical separation parameters as selected from the quality risk management evaluation, only the parameter tG and Flow have a significant and positive influence on critical resolution. The influence of T, Start %B, End %B or pH is small and probably not significant.

As a consequence the parameters tG and Start %B are confirmed to be critical separation parameters.


T	13.5 ± 1	[°C]	Start-%B	30 ± 1	%B
pH	2.1 ± 0.1		tG	8.1 ± 1	[min]
Flow Rate	0.5 ± 0.1	[mL/min]	End-%B	90 ± 1	%B
			Dwell Volume	0.12 ± 0	[mL]

Required Resolution	1.5
Successful Experiments	729
Success Rate	100 %
No of Factors	6
No of Levels	3
No of Experiments	729

Frequency Distribution




Regression Coefficients



Confirmation Run vs Predicted Run

Predicted Run



Confirmation Run





Table Confirmation vs Prediction

#	Name	Experimental tR [min]	Predicted tR [min]	Difference [min]	Error [%]
1	ImpD	0.84	1.10	-0.26	-23.95
2	ImpF	0.95	1.22	-0.27	-22.28
3	Amlodipine	1.23	1.51	-0.27	-18.09
4	ImpE	1.56	1.82	-0.26	-14.34
5	ImpB	2.44	2.69	-0.25	-9.15
6	ImpG	2.59	2.87	-0.28	-9.85
7	ImpH	2.72	2.96	-0.23	-7.82
8	ImpA	4.79	4.90	-0.11	-2.15

Chart Confirmation vs Prediction



Conclusions

In reference to the analytical target profile (ATP) the selected working point shows a separation off all impurities from each other and from the drug substance with more than resolution of 1.5 (2.6). The separation is good enough to guarantee the quantification limit 0.05%. The new proposed method has a run time minor than 6 min and it is compared to the old method (60) 10 times faster.

References and Comments

[1] Table 1 in Reference 1.

References:

1.- Exploring better column selectivity choices in ultra-high performance liquid chromatography using Quality by Design principles, R. Kormány, I. Molnár, H. J. Rieger /Journal of Pharmaceutical and Biomedical Analysis 80 (2013) 79-88

Model and document created with DryLab® version 4.3.1

Annex

Comparison

